

## Certificate of Analysis

**Product Category** 

**Blood Glucose Test Strips** 

Product Part Number

99GS700R51

Lot Number

212012002

Items	Specification	Analysis Data				
Precision testing 1						
(Test by Normal Control	CV≦5%	1.6%				
Solution)						
Precision testing 2						
(Test by venous blood /	CV≦5%	1.5%				
90±10 mg/dL)						
		1. $R^2 = 0.9995$				
	1. The correlation coefficient (R <sup>2</sup> )	2. Recovery ratio:				
Linearity testing	must be greater than 0.99.	1) 95.1% 2) 108.9%				
	2. The recovery of each sample	3) 109.7% 4) 113.9%				
	must be between 80%~120%.	5)_112.2%				

Final Product Evaluation: PASS [ V ] FAIL [	Final	Product	Evaluation:	PASS	[	٧	1	FAIL 【	1
---	-------	---------	-------------	------	---	---	---	--------	---

Storage Condition: Store at 4-30℃

Shelf-Life: 24 months

Complying Standards: ISO 13485: 2016 and 98/79/EC and ISO 15197:2013

Dates of Manufacturing:

2020.01.20

Date Issued:

2020.01.21

Expiration Date: 2022.01.18

**QA Supervisor** 

**QA Manager** 

**BIONIME CORPORATION** 

No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan

P 886 4 23692388 F 886 4 22617586

No.: FT-82-03-0024-0004F