

Certificate of Analysis

Product Category Blood Glucose Test Strips

Product Part Number 99GS700R21

Lot Number 2122B1301

Items	Specification	Analysis Data				
Precision testing 1 (Test by Normal Control	CV <u>≤</u> 5%	1.9%				
Solution)						
Precision testing 2						
(Test by venous blood /	CV <u>≦</u> 5%	2.0%				
90±10 mg/dL)						
Linearity testing	 The correlation coefficient (R²) must be greater than 0.99. The recovery of each sample must be between 80%~120%. 	1. R ² = 0.9973 2. Recovery ratio: 1) 93.4% 3) 100.6% 5) 110.9% 2) 97.2% 4) 101.6%				

Final	Product	Evalua	tion: P	ASS		٧	1	FAIL]
-------	---------	--------	---------	-----	--	---	---	------	--	---

Storage Condition: Store at 4-30°C

Shelf-Life: 24 months

Complying Standards: ISO 13485: 2016 and 98/79/EC and ISO 15197:2013

Dates of Manufacturing: 2022.11.14

Date Issued: 2022.11.15

Expiration Date: 2024.11.13

No.: FT-82-03-0024-0004F

QA Supervisor

QA Manager

No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan P 886 4 23692388 F 886 4 22617586