

## Certificate of Analysis

Product Category Blood Glucose Test Strips

Product Part Number 99GS770U11

Lot Number 2122A1211

ltems	Specification	Analysis Data 1.8%			
Precision testing 1 (Test by Normal Control	CV≦5%				
Solution)	,				
Precision testing 2					
(Test by venous blood /	CV <u>≤</u> 5%	1.7%			
90±10 mg/dL)					
Linearity testing	<ol> <li>The correlation coefficient (R²) must be greater than 0.99.</li> <li>The recovery of each sample must be between 80%~120%.</li> </ol>	1. R <sup>2</sup> = 1.0000 2. Recovery ratio: 1) 91.0% 2) 105.9% 3) 110.5% 4) 111.3% 5) 111.5%			

Final Product Evaluation: PASS		٧	1	FAIL 【	1
--------------------------------	--	---	---	--------	---

Storage Condition: Store at 4-30℃

Shelf-Life: 24 months

Complying Standards: ISO 13485: 2016 and 98/79/EC and ISO 15197:2013

Dates of Manufacturing: 2022.10.13

Date Issued: 2022.10.14

Expiration Date: 2024.10.12

No.: FT-82-03-0024-0004F

QA Supervisor

**QA** Manager

Carl Shill

**BIONIME CORPORATION** 

No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan P 886 4 23692388 F 886 4 22617586