

Certificate of Analysis

Product Category

Blood Glucose Test Strips

Product Part Number

99GS570U21

Lot Number

212330812

Items	Specification	Analysis Data
Precision testing 1 (Test by Normal Control Solution)	CV <u>≤</u> 5%	1.7%
Precision testing 2 (Test by venous blood / 90±10 mg/dL)	CV <u>≤</u> 5%	1.6%
Linearity testing	 The correlation coefficient (R²) must be greater than 0.99. The recovery of each sample must be between 80%~120%. 	1. $R^2 = 0.9979$ 2. Recovery ratio: 1) 103.7% 3) 114.4% 5) 111.7% 2) 106.9% 4) 116.0%

Final Product Evaluation: PASS [V] FAIL [1

Shelf-Life: 24 months

Complying Standards: ISO 13485: 2016 and 98/79/EC and ISO 15197:2013

Dates of Manufacturing:

2023.03.08

Date Issued:

2023.03.15

Expiration Date: 2025.03.07

QA Supervisor

QA Manager

BIONIME CORPORATION

No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan P 886 4 23692388 F 886 4 22617586

No.: FT-82-03-0024-0004F