

Certificate of Analysis

Product Category

Blood Glucose Test Strips

Product Part Number

99GS700R21

Lot Number

212331502

Items	Specification	Analysis Data			
Precision testing 1 (Test by Normal Control Solution)	CV <u>≤</u> 5%	1.6%			
Precision testing 2 (Test by venous blood / 90±10 mg/dL)	CV <u>≦</u> 5%	1.6%			
Linearity testing	 The correlation coefficient (R²) must be greater than 0.99. The recovery of each sample must be between 80%~120%. 	1. R ² = 0.9997 2. Recovery ratio: 1) 95.3% 3) 104.1% 5) 108.8% 2) 106.0% 4) 104.6%			

Final Product Evaluation: PASS		٧	1	FAIL 【	1
--------------------------------	--	---	---	--------	---

Storage Condition: Store at 4-30°C

Shelf-Life: 24 months

Complying Standards: ISO 13485: 2016 and 98/79/EC and ISO 15197:2013

Dates of Manufacturing:

2023.03.16

Date Issued:

2023.03.21

Expiration Date: 2025.03.15

QA Supervisor

QA Manager

Carl. Shih

BIONIME CORPORATION

No. 100, Sec. 2, Daqing St., South Dist., Taichung City 40242, Taiwan P 886 4 23692388 F 886 4 22617586