

Certificate of Analysis

Product Category Control Solution

Level Normal

Product Part Number 705-2GC550-050

Lot Number 2YD01A

Items	Specification	Analysis Data
Accuracy testing (Glucose level test by YSI 2300)	96 mg/dL $\pm$ 4 mg/dL	Test mean : 95.8 mg/dL
Within Lot Precision (CV %)	CV % $\leq$ 4%	CV % = 0.7 %

Final Product Evaluation : PASS 【 **V** 】 FAIL 【 】

Storage Condition : Store at 2-30°C

Shelf-Life : 24 months

Complying Standards : ISO 13485:2016 and 98/79/EC and ISO 15197:2013

Dates of Manufacturing: 2025.02.05

Date Issued: 2025.04.29

Expiration Date: 2027.02.04

Guo Gu
QA Supervisor

Carl Shih
QA Manager

No.:FT-82-03-0024-0007E

BIONIME CORPORATION
No. 100, Sec. 2, Daqing St., South Dist.,
Taichung City 40242, Taiwan
P 886 4 23692388 F 886 4 22617586